

# Chapter 1: Introduction



FIGURE 1-1 Model 335 front view

## 1.1 Product Description

### Features:

- Operates down to 300 mK with appropriate NTC RTD sensors
- Two sensor inputs
- Two configurable PID control loops providing 50 W and 25 W or 75 W and 1 W
- Autotuning automatically calculates PID parameters
- Automatically switch sensor inputs using zones to allow continuous measurement and control from 300 mK to 1505 K
- Custom display setup allows you to label each sensor input
- USB and IEEE-488 interfaces
- Supports diode, RTD, and thermocouple temperature sensors
- Sensor excitation current reversal eliminates thermal EMF errors for resistance sensors
- $\pm 10$  V analog voltage outputs, alarms, and relays

Designed with the user and ease of use in mind, the Model 335 temperature controller offers many user-configurable features and advanced functions that until now have been reserved for more expensive, high-end temperature controllers. The Model 335 is the first two-channel temperature controller available with user configurable heater outputs delivering a total of 75 W of low noise heater power: 50 W and 25 W, or 75 W and 1 W. With that much heater power packed into an affordable half-rack sized instrument, the Model 335 gives you more power and control than ever.

Control outputs are equipped with both hardware and software features allowing you, and not your temperature controller, to easily control your experiments. Output one functions as a current output while output two can be configured in either current or voltage mode. With output two in voltage mode, it functions as a  $\pm 10$  V analog output while still providing 1 W of heater power and full closed loop proportional-integral-derivative (PID) control capability. Alarms and relays are included to help automate secondary control functions. The improved autotuning feature of the Model 335 can be used to automatically calculate PID control parameters, so you spend less time tuning your controller and more time conducting experiments.

The Model 335 supports the industry's most advanced line of cryogenic temperature sensors as manufactured by Lake Shore, including diodes, resistance temperature detectors (RTDs), and thermocouples. The controller's zone tuning feature allows you to measure and control temperatures seamlessly from 300 mK to over 1,500 K. This feature automatically switches temperature sensor inputs when your temperature range goes beyond the useable range of a given sensor. You'll never again have to be concerned with temperature sensor over or under errors and measurement continuity issues.

The intuitive front panel layout and keypad logic, bright vacuum fluorescent display, and LED indicators enhance the user-friendly front panel interface of the Model 335. Four standard display modes are offered to accommodate different instrument configurations and user preferences. Say goodbye to sticky notes and hand written labels, as the ability to custom label sensor inputs eliminates the guesswork in remembering or determining the location to which a sensor input is associated. These features, combined with USB and IEEE-488 interfaces and intuitive menu structure and logic supports efficiency and ease of use.

As a replacement to our popular Model 331 and 332 temperature controllers, the Model 335 offers software emulation modes for literal drop-in compatibility. The commands you are accustomed to sending to the Model 331 and 332 will either be interpreted directly or translated to the most appropriate Model 335 setting. The Model 335 comes standard-equipped with all of the functionality of the controllers it replaces, but offers additional features that save you time and money. With the Model 335, you get a temperature controller you control from the world leader in cryogenic thermometry.

### 1.1.1 Sensor Inputs

The Model 335 offers two standard sensor inputs that are compatible with diode and RTD temperature sensors. The field-installable Model 3060 option adds thermocouple functionality to both inputs.

Sensor inputs feature a high-resolution 24-bit analog-to-digital converter and each of the two powered outputs function as separate current sources. Both sensor inputs are optically isolated from other circuits to reduce noise and to deliver repeatable sensor measurements. Current reversal eliminates thermal electromagnetic field (EMF) errors in resistance sensors. Ten excitation currents facilitate temperature measurement and control down to 300 mK using appropriate negative temperature coefficient (NTC) RTDs. Autorange mode automatically scales excitation current in NTC RTDs to reduce self heating at low temperatures as sensor resistance changes by many orders of magnitude. Temperatures down to 1.4 K can be measured and controlled using silicon or GaAlAs diodes. Software selects the appropriate excitation current and signal gain levels when the sensor type is entered via the instrument front panel. To increase your productivity, the unique zone setting feature automatically switches sensor inputs, enabling you to measure temperatures from 300 mK to over 1,500 K without interrupting your experiment.

The Model 335 includes standard temperature sensor response curves for silicon diodes, platinum RTDs, ruthenium oxide RTDs, and thermocouples. Non-volatile memory can also store up to 39 200-point CalCurves for Lake Shore calibrated temperature sensors or user curves. A built-in SoftCal algorithm can be used to generate curves for silicon diodes and platinum RTDs that can be stored as user curves. Temperature sensor calibration data can be easily loaded into the Model 335 temperature controller and manipulated using the Lake Shore Curve Handler™ software program.

## 1.1.2 Temperature Control

Providing a total of 75 W of heater power, the Model 335 is the most powerful half rack temperature controller available. Designed to deliver very clean heater power, precise temperature control is ensured throughout your full scale temperature range for excellent measurement reliability, efficiency and throughput. Two independent PID control outputs can be configured to supply 50 W and 25 W or 75 W and 1 W of heater power. Precise control output is calculated based on your temperature set-point and feedback from the control sensor. Wide tuning parameters accommodate most cryogenic cooling systems and many high-temperature ovens commonly used in laboratories. PID values can be manually set for fine control or the improved auto-tuning feature can automate the tuning process.

The Model 335 autotuning method calculates PID parameters and provides feedback to help build zone tables. The setpoint ramp feature provides smooth, continuous set-point changes and predictable approaches to setpoint without the worry of overshoot or excessive settling times. The instrument's zone tuning feature automatically switches temperature sensor inputs when your temperature range goes beyond the useable range of a given sensor. This feature combined with the instrument's ability to scale the sensor excitation through ten pre-loaded current settings allows the Model 335 to provide continuous measurement and control from 300 mK to 1505 K.

Both control outputs are variable DC current sources referenced to chassis ground. As a factory default, Outputs 1 and 2 provide 50 W and 25 W of continuous power respectively, both to a 50  $\Omega$  or 25  $\Omega$  load. For increased functionality, Output 2 can also be set to voltage mode. When set to voltage mode, it functions as a  $\pm 10$  V analog output while still providing 1 W of heater power and full closed loop PID control capability. While in this mode, output 1 can provide up to 75 W of heater power to a 25  $\Omega$  load.

Temperature limit settings for inputs are provided as a safeguard against system damage. Each input is assigned a temperature limit, and if any input exceeds that limit, both control channels are automatically disabled.

## 1.1.3 Interface

The Model 335 is standard equipped with universal serial bus (USB) and parallel (IEEE-488) interfaces. In addition to gathering data, nearly every function of the instrument can be controlled via computer interface. You can download the Lake Shore Curve Handler™ software program to your computer to easily enter and manipulate sensor calibration curves for storage in the instrument's non-volatile memory.

The USB interface emulates an RS-232C serial port at a fixed 57,600 baud rate, but with the physical plug-ins of a USB. It also allows you to download firmware upgrades, ensuring the most current firmware version is loaded into your instrument without having to physically change your instrument.

Both sensor inputs are equipped with a high and low alarm which offers latching and non-latching operation. The two relays can be used in conjunction with the alarms to alert you of a fault condition and perform simple on-off control. Relays can be assigned to any alarm or operated manually.

The  $\pm 10$  V analog voltage output can be configured to send a voltage proportional to temperature to a strip chart recorder or data acquisition system. You may select the scale and data sent to the output, including temperature or sensor units.




FIGURE 1-2 Model 335 rear panel

### 1.1.4 Configurable Display

The Model 335 offers a bright, vacuum fluorescent display that simultaneously displays up to four readings. You can display both control loops, or if you need to monitor just one input, you can display just that one in greater detail. Or you can custom configure each display location to suit your experiment. Data from any input can be assigned to any of the locations, and your choice of temperature sensor units can be displayed. For added convenience, you can also custom label each sensor input, eliminating the guesswork in remembering or determining the location to which a sensor input is associated.

- Two Input/Output Display with Labels: Standard display option featuring two inputs and associated outputs.
- Custom Display with Labels: Reading locations can be user configured to accommodate application needs. Here, the input names are shown above the measurement readings along with the designated input letters.
- Intuitive Menu Structure: Logical navigation allows you to spend more time on research and less time on setup.



FIGURE 1-3 Displays showing two input/output display with labels, custom display with labels and the intuitive menu structure

### 1.1.5 Model 3060 Thermocouple Input Option

The field installable Model 3060 thermocouple input option adds thermocouple functionality to both inputs. While the option can be easily removed, this is not necessary as the standard inputs remain fully functional when they are not being used to measure thermocouple temperature sensors. Calibration for the option is stored on the card so it can be installed in the field and used with multiple Model 335 temperature controllers without recalibration.

## 1.2 Sensor Selection

Silicon diodes are the best choice for general cryogenic use from 1.4 K to above room temperature. Diodes are economical to use because they follow a standard curve and are interchangeable in many applications. They are not suitable for use in ionizing radiation or magnetic fields.

Cernox™ thin-film RTDs offer high sensitivity and low magnetic field-induced errors over the 0.3 K to 420 K temperature range. Cernox sensors require calibration.

Platinum RTDs offer high uniform sensitivity from 30 K to over 800 K. With excellent reproducibility, they are useful as thermometry standards. They follow a standard curve above 70 K and are interchangeable in many applications.

		Model	Useful range	Magnetic field use
<b>Diodes</b>	Silicon Diode	DT-670-SD	1.4 K to 500 K	$T \geq 60 \text{ K} \ \& \ B \leq 3 \text{ T}$
	Silicon Diode	DT-670E-BR	30 K to 500 K	$T \geq 60 \text{ K} \ \& \ B \leq 3 \text{ T}$
	Silicon Diode	DT-414	1.4 K to 375 K	$T \geq 60 \text{ K} \ \& \ B \leq 3 \text{ T}$
	Silicon Diode	DT-421	1.4 K to 325 K	$T \geq 60 \text{ K} \ \& \ B \leq 3 \text{ T}$
	Silicon Diode	DT-470-SD	1.4 K to 500 K	$T \geq 60 \text{ K} \ \& \ B \leq 3 \text{ T}$
	Silicon Diode	DT-471-SD	10 K to 500 K	$T \geq 60 \text{ K} \ \& \ B \leq 3 \text{ T}$
	GaAlAs Diode	TG-120-P	1.4 K to 325 K	$T > 4.2 \text{ K} \ \& \ B \leq 5 \text{ T}$
	GaAlAs Diode	TG-120-PL	1.4 K to 325 K	$T > 4.2 \text{ K} \ \& \ B \leq 5 \text{ T}$
	GaAlAs Diode	TG-120-SD	1.4 K to 500 K	$T > 4.2 \text{ K} \ \& \ B \leq 5 \text{ T}$
<b>Positive Temperature Coefficient RTDs</b>	100 $\Omega$ Platinum	PT-102/3	14 K to 873 K	$T > 40 \text{ K} \ \& \ B \leq 2.5 \text{ T}$
	100 $\Omega$ Platinum	PT-111	14 K to 673 K	$T > 40 \text{ K} \ \& \ B \leq 2.5 \text{ T}$
	Rhodium-Iron	RF-800-4	1.4 K to 500 K	$T > 77 \text{ K} \ \& \ B \leq 8 \text{ T}$
	Rhodium-Iron	RF-100T/U	1.4 K to 325 K	$T > 77 \text{ K} \ \& \ B \leq 8 \text{ T}$
<b>Negative Temperature Coefficient RTDs</b>	Cernox™	CX-1010	0.3 K to 325 K <sup>1</sup>	$T > 2 \text{ K} \ \& \ B \leq 19 \text{ T}$
	Cernox™	CX-1030-HT	0.3 K to 420 K <sup>1,3</sup>	$T > 2 \text{ K} \ \& \ B \leq 19 \text{ T}$
	Cernox™	CX-1050-HT	1.4 K to 420 K <sup>1</sup>	$T > 2 \text{ K} \ \& \ B \leq 19 \text{ T}$
	Cernox™	CX-1070-HT	4 K to 420 K <sup>1</sup>	$T > 2 \text{ K} \ \& \ B \leq 19 \text{ T}$
	Cernox™	CX-1080-HT	20 K to 420 K <sup>1</sup>	$T > 2 \text{ K} \ \& \ B \leq 19 \text{ T}$
	Germanium	GR-300-AA	0.35 K to 100 K <sup>3</sup>	Not recommended
	Germanium	GR-1400-AA	1.8 K to 100 K <sup>3</sup>	Not recommended
	Carbon-Glass	CGR-1-500	1.4 K to 325 K	$T > 2 \text{ K} \ \& \ B \leq 19 \text{ T}$
	Carbon-Glass	CGR-1-1000	1.7 K to 325 K <sup>2</sup>	$T > 2 \text{ K} \ \& \ B \leq 19 \text{ T}$
	Carbon-Glass	CGR-1-2000	2 K to 325 K <sup>2</sup>	$T > 2 \text{ K} \ \& \ B \leq 19 \text{ T}$
	Rox™	RX-102	0.3 K to 40 K <sup>3</sup>	$T > 2 \text{ K} \ \& \ B \leq 10 \text{ T}$
	Rox™	RX-103	1.4 K to 40 K	$T > 2 \text{ K} \ \& \ B \leq 10 \text{ T}$
	Rox™	RX-202	0.3 K to 40 K <sup>3</sup>	$T > 2 \text{ K} \ \& \ B \leq 10 \text{ T}$
<b>Thermocouples 3060</b>	Type K	9006-006	3.2 K to 1505 K	Not recommended
	Type E	9006-004	3.2 K to 934 K	Not recommended
	Chromel-AuFe 0.07%	9006-002	1.2 K to 610 K	Not recommended

<sup>1</sup> Non-HT version maximum temperature: 325 K

<sup>2</sup> Low temperature limited by input resistance range

<sup>3</sup> Low temperature specified with self-heating error:  $\leq 5 \text{ mK}$

TABLE 1-1 Sensor temperature range

	Example Lake Shore Sensor	Temperature	Nominal Resistance/Voltage	Typical Sensor Sensitivity <sup>4</sup>	Measurement Resolution: Temperature Equivalents	Electronic Accuracy: Temperature Equivalents	Temperature Accuracy including Electronic Accuracy, CalCurve and Calibrated Sensor	Electronic Control Stability <sup>5</sup> : Temperature Equivalents
Silicon Diode	DT-670-CO-13 with 1.4H calibration	1.4 K	1.664 V	-12.49 mV/K	0.8 mK	±13 mK	±25 mK	±1.6 mK
		77 K	1.028 V	-1.73 mV/K	5.8 mK	±76 mK	±98 mK	±11.6 mK
		300 K	0.5597 V	-2.3 mV/K	4.4 mK	±47 mK	±79 mK	±8.8 mK
		500 K	0.0907 V	-2.12 mV/K	4.7 mK	±40 mK	±90 mK	±9.4 mK
Silicon Diode	DT-470-SD-13 with 1.4H calibration	1.4 K	1.6981 V	-13.1 mV/K	0.8 mK	±13 mK	±25 mK	±1.6 mK
		77 K	1.0203 V	-1.92 mV/K	5.2 mK	±69 mK	±91 mK	±10.4 mK
		300 K	0.5189 V	-2.4 mV/K	4.2 mK	±45 mK	±77 mK	±8.4 mK
		475 K	0.0906 V	-2.22 mV/K	4.5 mK	±38 mK	±88 mK	±9 mK
GaAlAs Diode	TG-120-SD with 1.4H calibration	1.4 K	5.391 V	-97.5 mV/K	0.2 mK	±7 mK	±19 mK	±0.4 mK
		77 K	1.422 V	-1.24 mV/K	16 mK	±180 mK	±202 mK	±32 mK
		300 K	0.8978 V	-2.85 mV/K	7 mK	±60 mK	±92 mK	±14 mK
		475 K	0.3778 V	-3.15 mV/K	6.4 mK	±38 mK	±88 mK	±13 mK
100 Ω Platinum RTD 500 Ω Full Scale	PT-103 with 14J calibration	30 K	3.660 Ω	0.191 Ω/K	1.1 mK	±13 mK	±23 mK	±2.2 mK
		77 K	20.38 Ω	0.423 Ω/K	0.5 mK	±10 mK	±22 mK	±1.0 mK
		300 K	110.35 Ω	0.387 Ω/K	5.2 mK	±39 mK	±62 mK	±10.4 mK
		500 K	185.668 Ω	0.378 Ω/K	5.3 mK	±60 mK	±106 mK	±10.6 mK
Cernox™	CX-1010-SD with 0.3L calibration	0.3 K	2322.4 Ω	-10785 Ω/K	8.5 μK	±0.1 mK	±3.6 mK	±17 μK
		0.5 K	1248.2 Ω	-2665.2 Ω/K	26 μK	±0.2 mK	±4.7 mK	±52 μK
		4.2 K	277.32 Ω	-32.209 Ω/K	140 μK	±3.8 mK	±8.8 mK	±280 μK
		300 K	30.392 Ω	-0.0654 Ω/K	23 mK	±339 mK	±414 mK	±46 mK
Cernox™	CX-1050-SD-HT <sup>6</sup> with 1.4M calibration	1.4 K	26566 Ω	-48449 Ω/K	20 μK	±0.3 mK	±5.3 mK	±40 μK
		4.2 K	3507.2 Ω	-1120.8 Ω/K	196 μK	±2.1 mK	±7.1 mK	±392 μK
		77 K	205.67 Ω	-2.4116 Ω/K	1.9 mK	±38 mK	±54 mK	±3.8 mK
		420 K	45.03 Ω	-0.0829 Ω/K	18 mK	±338 mK	±403 mK	±36 mK
Germanium	GR-300-AA with 0.3D calibration	0.35 K	18225 Ω	-193453 Ω/K	4 μK	±48 μK	±4.2 mK	±8 μK
		1.4 K	449 Ω	-581 Ω/K	41 μK	±481 μK	±4.7 mK	±82 μK
		4.2 K	94 Ω	-26.6 Ω/K	56 μK	±1.8 mK	±6.8 mK	±112 μK
		100 K	2.7 Ω	-0.024 Ω/K	6.3 mK	±152 mK	±175 mK	±12.6 mK
Germanium	GR-1400-AA with 1.4D calibration	1.8 K	15288 Ω	-26868 Ω/K	28 μK	±302 μK	±4.5 mK	±56 μK
		4.2 K	1689 Ω	-862 Ω/K	91 μK	±900 μK	±5.1 mK	±182 μK
		10 K	253 Ω	-62.0 Ω/K	73 μK	±1.8 mK	±6.8 mK	±146 μK
		100 K	2.8 Ω	-0.021 Ω/K	7.1 mK	±177 mK	±200 mK	±14.2 mK
Carbon-Glass	CGR-1-500 with 1.4L calibration	1.4 K	103900 Ω	-520000 Ω/K	13 μK	±0.1 mK	±4.1 mK	±26 μK
		4.2 K	584.6 Ω	-422.3 Ω/K	63 μK	±0.8 mK	±4.8 mK	±126 μK
		77 K	14.33 Ω	-0.098 Ω/K	4.6 mK	±108 mK	±133 mK	±9.2 mK
		300 K	8.55 Ω	-0.0094 Ω/K	16 mK	±760 mK	±865 mK	±32 mK
Rox™	RX-102A-AA with 0.3B calibration	0.5 K	3701 Ω	-5478 Ω/K	41 μK	±0.5 mK	±5 mK	±82 μK
		1.4 K	2005 Ω	-667 Ω/K	128 μK	±1.4 mK	±6.4 mK	±256 μK
		4.2 K	1370 Ω	-80.3 Ω/K	902 μK	±8 mK	±24 mK	±1.8 mK
		40 K	1049 Ω	-1.06 Ω/K	62 mK	±500 mK	±537 K	±124 mK
Thermocouple 50 mV 3060-F	Type K	75 K	-5862.9 μV	15.6 μV/K	26 mK	±0.25 K <sup>7</sup>	Calibration not available from Lake Shore	±52 mK
		300 K	1075.3 μV	40.6 μV/K	10 mK	±0.038 K <sup>7</sup>		±20 mK
		600 K	13325 μV	41.7 μV/K	10 mK	±0.184 K <sup>7</sup>		±20 mK
		1505 K	49998.3 μV	36.006 μV/K	11 mK	±0.73 K <sup>7</sup>		±22 mK

<sup>4</sup> Typical sensor sensitivities were taken from representative calibrations for the sensor listed

<sup>5</sup> Control stability of the electronics only, in an ideal thermal system

<sup>6</sup> Non-HT version maximum temperature: 325 K

<sup>7</sup> Accuracy specification does not include errors from room temperature compensation

TABLE 1-2 Typical sensor performance

## 1.3 Model 335 Specifications

Full specifications are provided on our website. Please see:  
<https://www.lakeshore.com/335>.

### 1.3.1 Input Specifications

	Sensor temperature coefficient	Input range	Excitation current	Display resolution	Measurement resolution	Electronic accuracy (at 25 °C)	Measurement temperature coefficient	Electronic stability <sup>1</sup>
<b>Diode</b>	Negative	0 V to 2.5 V	10 $\mu$ A $\pm 0.05\%$ <sup>2,3</sup>	100 $\mu$ V	10 $\mu$ V	$\pm 80 \mu\text{V} \pm 0.005\%$ of rdg	(10 $\mu$ V + 0.0005% of rdg)/°C	$\pm 20 \mu\text{V}$
		0 V to 10 V	10 $\mu$ A $\pm 0.05\%$ <sup>2,3</sup>	1 mV	20 $\mu$ V	$\pm 320 \mu\text{V} \pm 0.01\%$ of rdg	(20 $\mu$ V + 0.0005% of rdg)/°C	$\pm 40 \mu\text{V}$
<b>PTC RTD</b>	Positive	0 $\Omega$ to 10 $\Omega$	1 mA <sup>4</sup>	1 m $\Omega$	0.2 m $\Omega$	$\pm 0.002 \Omega \pm 0.01\%$ of rdg	(0.01 m $\Omega$ + 0.001% of rdg)/°C	$\pm 0.4 \text{ m}\Omega$
		0 $\Omega$ to 30 $\Omega$	1 mA <sup>4</sup>	1 m $\Omega$	0.2 m $\Omega$	$\pm 0.002 \Omega \pm 0.01\%$ of rdg	(0.03 m $\Omega$ + 0.001% of rdg)/°C	$\pm 0.4 \text{ m}\Omega$
		0 $\Omega$ to 100 $\Omega$	1 mA <sup>4</sup>	10 m $\Omega$	2 m $\Omega$	$\pm 0.004 \Omega \pm 0.01\%$ of rdg	(0.1 m $\Omega$ + 0.001% of rdg)/°C	$\pm 4 \text{ m}\Omega$
		0 $\Omega$ to 300 $\Omega$	1 mA <sup>4</sup>	10 m $\Omega$	2 m $\Omega$	$\pm 0.004 \Omega \pm 0.01\%$ of rdg	(0.3 m $\Omega$ + 0.001% of rdg)/°C	$\pm 4 \text{ m}\Omega$
		0 $\Omega$ to 1 k $\Omega$	1 mA <sup>4</sup>	100 m $\Omega$	20 m $\Omega$	$\pm 0.04 \Omega \pm 0.02\%$ of rdg	(1 m $\Omega$ + 0.001% of rdg)/°C	$\pm 40 \text{ m}\Omega$
		0 $\Omega$ to 3 k $\Omega$	1 mA <sup>4</sup>	100 m $\Omega$	20 m $\Omega$	$\pm 0.04 \Omega \pm 0.02\%$ of rdg	(3 m $\Omega$ + 0.001% of rdg)/°C	$\pm 40 \text{ m}\Omega$
		0 $\Omega$ to 10 k $\Omega$	1 mA <sup>4</sup>	1 $\Omega$	200 m $\Omega$	$\pm 0.4 \Omega \pm 0.02\%$ of rdg	(10 m $\Omega$ + 0.001% of rdg)/°C	$\pm 400 \text{ m}\Omega$
<b>NTC RTD 10 mV</b>	Negative	0 $\Omega$ to 10 $\Omega$	1 mA <sup>4</sup>	1 m $\Omega$	0.15 m $\Omega$	$\pm 0.002 \Omega \pm 0.06\%$ of rdg	(0.01 m $\Omega$ + 0.001% of rdg)/°C	$\pm 0.3 \text{ m}\Omega$
		0 $\Omega$ to 30 $\Omega$	300 $\mu$ A <sup>4</sup>	1 m $\Omega$	0.45 m $\Omega$	$\pm 0.002 \Omega \pm 0.06\%$ of rdg	(0.03 m $\Omega$ + 0.0015% of rdg)/°C	$\pm 0.9 \text{ m}\Omega$
		0 $\Omega$ to 100 $\Omega$	100 $\mu$ A <sup>4</sup>	10 m $\Omega$	1.5 m $\Omega$	$\pm 0.01 \Omega \pm 0.04\%$ of rdg	(0.1 m $\Omega$ + 0.001% of rdg)/°C	$\pm 3 \text{ m}\Omega$
		0 $\Omega$ to 300 $\Omega$	30 $\mu$ A <sup>4</sup>	10 m $\Omega$	4.5 m $\Omega$	$\pm 0.01 \Omega \pm 0.04\%$ of rdg	(0.3 m $\Omega$ + 0.0015% of rdg)/°C	$\pm 9 \text{ m}\Omega$
		0 $\Omega$ to 1 k $\Omega$	10 $\mu$ A <sup>4</sup>	100 m $\Omega$	15 m $\Omega$ + 0.002% of rdg	$\pm 0.1 \Omega \pm 0.04\%$ of rdg	(1 m $\Omega$ + 0.001% of rdg)/°C	$\pm 30 \text{ m}\Omega \pm 0.004\%$ of rdg
		0 $\Omega$ to 3 k $\Omega$	3 $\mu$ A <sup>4</sup>	100 m $\Omega$	45 m $\Omega$ + 0.002% of rdg	$\pm 0.1 \Omega \pm 0.04\%$ of rdg	(3 m $\Omega$ + 0.0015% of rdg)/°C	$\pm 90 \text{ m}\Omega \pm 0.004\%$ of rdg
		0 $\Omega$ to 10 k $\Omega$	1 $\mu$ A <sup>4</sup>	1 $\Omega$	150 m $\Omega$ + 0.002% of rdg	$\pm 1.0 \Omega \pm 0.04\%$ of rdg	(10 m $\Omega$ + 0.001% of rdg)/°C	$\pm 300 \text{ m}\Omega \pm 0.004\%$ of rdg
		0 $\Omega$ to 30 k $\Omega$	300 nA <sup>4</sup>	1 $\Omega$	450 m $\Omega$ + 0.002% of rdg	$\pm 2.0 \Omega \pm 0.04\%$ of rdg	(30 m $\Omega$ + 0.0015% of rdg)/°C	$\pm 900 \text{ m}\Omega \pm 0.004\%$ of rdg
		0 $\Omega$ to 100 k $\Omega$	100 nA <sup>4</sup>	10 $\Omega$	1.5 $\Omega$ + 0.005% of rdg	$\pm 10.0 \Omega \pm 0.04\%$ of rdg	(100 m $\Omega$ + 0.002% of rdg)/°C	$\pm 3 \Omega \pm 0.01\%$ of rdg
<b>Thermocouple</b>	Positive	$\pm 50 \text{ mV}$	NA	1 $\mu$ V	0.4 $\mu$ V	$\pm 1 \mu\text{V} \pm 0.05\%$ of rdg <sup>5</sup>	(0.1 $\mu$ V + 0.001% of rdg)/°C	$\pm 0.8 \mu\text{V}$

<sup>1</sup> Control stability of the electronics only, in ideal thermal system

<sup>2</sup> Current source error has negligible effect on measurement accuracy

<sup>3</sup> Diode input excitation can be set to 1 mA

<sup>4</sup> Current source error is removed during calibration

<sup>5</sup> Accuracy specification does not include errors from room temperature compensation

TABLE 1-3 Input specifications

### 1.3.2 Sensor Input Configuration

	Diode/RTD	Thermocouple
<b>Measurement type</b>	4-lead differential	2-lead differential, room temperature compensated
<b>Excitation</b>	Constant current with current reversal for RTDs	NA
<b>Supported sensors</b>	<b>Diodes:</b> Silicon, GaAlAs <b>RTDs:</b> 100 $\Omega$ Platinum, 1000 $\Omega$ Platinum Germanium, Carbon-Glass, Cernox™, and Rox™	Most thermocouple types
<b>Standard curves</b>	DT-470, DT-670, DT-500-D, DT-500-E1, PT-100, PT-1000, RX-102A, RX-202A	Type E, Type K, Type T, AuFe 0.07% vs. Cr, AuFe 0.03% vs. CR
<b>Input connector</b>	6-pin DIN	Screw terminals in a ceramic isothermal block

TABLE 1-4 *Sensor input configuration*

### 1.3.3 Thermometry

<b>Number of inputs</b>	2
<b>Input configuration</b>	Inputs can be configured from the front panel to accept any of the supported input types. Thermocouple inputs require an optional input card that can be installed in the field. Once installed, the thermocouple input can be selected from the front panel like any other input type.
<b>Isolation</b>	Sensor inputs optically isolated from other circuits but not each other
<b>A/D resolution</b>	24-bit
<b>Input accuracy</b>	Sensor dependent, refer to Input Specifications table
<b>Measurement resolution</b>	Sensor dependent, refer to Input Specifications table
<b>Maximum update rate</b>	10 rdg/s on each input, 5 rdg/s when configured as 100 k $\Omega$ NTC RTD with reversal on
<b>Autorange</b>	Automatically selects appropriate NTC RTD or PTC RTD range
<b>User curves</b>	Room for 39 200-point CalCurves™ or user curves
<b>SoftCal™</b>	Improves accuracy of DT-470 diode to $\pm 0.25$ K from 30 K to 375 K; improves accuracy of platinum RTDs to $\pm 0.25$ K from 70 K to 325 K; stored as user curves
<b>Math</b>	Maximum and minimum
<b>Filter</b>	Averages 2 to 64 input readings

### 1.3.4 Control

There are two control outputs.

#### 1.3.4.1 Heater Outputs (Outputs 1 and 2)

<b>Control type</b>	Closed loop digital PID with manual heater output or open loop; warm up mode (Output 2 only)
<b>Update rate</b>	10/s
<b>Tuning</b>	Autotune (one loop at a time), PID, PID zones
<b>Control stability</b>	Sensor dependent, see Input Specifications table
<b>PID control settings</b>	
<b>Proportional (gain)</b>	0 to 1000 with 0.1 setting resolution
<b>Integral (reset)</b>	1 to 1000 (1000/s) with 0.1 setting resolution
<b>Derivative (rate)</b>	1 to 200% with 1% resolution
<b>Manual output</b>	0 to 100% with 0.01% setting resolution
<b>Zone control</b>	10 temperature zones with P, I, D, manual heater out, heater range, control channel, ramp rate
<b>Setpoint ramping</b>	0.1 K/min to 100 K/min



Type	Variable DC current source		
Control modes	Closed loop digital PID with manual output or open loop		
D/A resolution	16-bit		
	25 $\Omega$ setting		50 $\Omega$ setting
Max power	75 W*	50 W	50 W
Max current	1.73 A	1.41 A	1 A
Voltage compliance (min)	43.3 V	35.4 V	50 V
Heater load for max power	25 $\Omega$	25 $\Omega$	50 $\Omega$
Heater load range	10 $\Omega$ to 100 $\Omega$		
Ranges	3 (decade steps in power)		
Heater noise	0.12 $\mu$ A RMS (dominated by line frequency and its harmonics)		
Heater connector	Dual banana		
Grounding	Output referenced to chassis ground		
Safety limits	Curve temperature, power up heater off, short circuit protection		

\* 75 W only available when Output 2 is in voltage mode

TABLE 1-5 Output 1

Type	Variable DC current source or voltage source		
	Current mode		Voltage mode
Control modes	Closed loop digital PID with manual output, zone, open loop		Closed loop digital PID with manual output, zone, open loop, warm up, monitor out
D/A resolution	15-bit		16-bit (bipolar)/15-bit (unipolar)
	25 $\Omega$ setting	50 $\Omega$ setting	N/A
Max power	25 W	25 W	1 W
Max current	1 A	0.71 A	100 mA
Voltage compliance (min)	25 V	35.4 V	$\pm 10$ V
Heater load for max power	25 $\Omega$	50 $\Omega$	100 $\Omega$
Heater load range	10 $\Omega$ to 100 $\Omega$		100 $\Omega$ min (short circuit protected)
Ranges	3 (decade steps in power)		N/A
Heater noise	0.12 $\mu$ A RMS		0.3 mV RMS
Heater connector	Dual banana		Detachable terminal block
Grounding	Output referenced to chassis ground		
Safety limits	Curve temperature, power up heater off, short circuit protection		

TABLE 1-6 Output 2

#### Warm up heater mode settings (Output 2 only)

Warm up percentage	0 to 100% with 1% resolution
Warm up mode	Continuous control or auto-off

#### Monitor Output settings (Output 2 voltage only)

Scale	User selected
Data source	Temperature or sensor units
Settings	Input, source, top of scale, bottom of scale or manual
Update rate	10/s
Range	$\pm 10$ V
Resolution	16-bit, 0.3 mV
Accuracy	$\pm 2.5$ mV
Noise	0.3 mV RMS
Minimum load resistance	100 $\Omega$ (short-circuit protected)
Connector	Detachable terminal block

### 1.3.5 Front Panel

<b>Display</b>	2-line by 20-character, 9 mm character height, vacuum fluorescent display
<b>Number of reading displays</b>	1 to 4
<b>Display units</b>	K, °C, V, mV, $\Omega$
<b>Reading source</b>	Temperature, sensor units, max, and min
<b>Display update rate</b>	2 rdg/s
<b>Temperature display resolution</b>	0.001° from 0° to 99.999°, 0.01° from 100° to 999.99°, 0.1° above 1000°
<b>Sensor units display resolution</b>	Sensor dependent, to 5 digits
<b>Other displays</b>	Sensor name, setpoint, heater range, heater output, and PID
<b>Setpoint setting resolution</b>	Same as display resolution (actual resolution is sensor dependent)
<b>Heater output display</b>	Numeric display in percent of full scale for power or current
<b>Heater output resolution</b>	1%
<b>Display annunciators</b>	Control input, alarm, tuning
<b>LED annunciators</b>	Remote, alarm, control outputs
<b>Keypad</b>	25-key silicone elastomer keypad
<b>Front panel features</b>	Front panel curve entry, display brightness control, and keypad lock-out

### 1.3.6 Interface

<b>IEEE-488.2</b>	
<b>Capabilities</b>	SH1, AH1, T5, L4, SR1, RL1, PP0, DC1, DT0, CO, E1
<b>Reading rate</b>	To 10 rdg/s on each input
<b>Software support</b>	LabVIEW™ driver (contact Lake Shore for availability)
<b>USB</b>	
<b>Function</b>	Emulates a standard RS-232 serial port
<b>Baud Rate</b>	57,600
<b>Connector</b>	B-type USB connector
<b>Reading rate</b>	To 10 rdg/s on each input
<b>Software support</b>	LabVIEW™ driver (contact Lake Shore for availability)
<b>Special interface features</b>	Model 331/332 command emulation mode
<b>Alarms</b>	
<b>Number</b>	2, high and low for each input
<b>Data source</b>	Temperature or sensor units
<b>Settings</b>	Source, high setpoint, low setpoint, deadband, latching or non-latching, audible on/off, and visible on/off
<b>Actuators</b>	Display annunciator, beeper, and relays
<b>Relays</b>	
<b>Number</b>	2
<b>Contacts</b>	Normally open (NO), normally closed (NC), and common (C)
<b>Contact rating</b>	30 VDC at 3 A
<b>Operation</b>	Activate relays on high, low, or both alarms for any input, or manual mode
<b>Connector</b>	Detachable terminal block

### 1.3.7 General

<b>Ambient temperature</b>	15 °C to 35 °C at rated specifications; 5 °C to 40 °C at reduced specifications
<b>Power requirement</b>	100, 120, 220, 240, VAC, $\pm 10\%$ , 50 or 60 Hz, 210 VA
<b>Size</b>	217 mm W $\times$ 90 mm H $\times$ 317 mm D (8.5 in $\times$ 3.5 in $\times$ 14.5 in), half rack
<b>Weight</b>	7.6 kg (16.8 lb)
<b>Approval</b>	CE mark

## 1.4 Safety Summary and Symbols

Observe these general safety precautions during all phases of instrument operation, service, and repair. Failure to comply with these precautions or with specific warnings elsewhere in this manual violates safety standards of design, manufacture, and intended instrument use. Lake Shore Cryotronics, Inc. assumes no liability for Customer failure to comply with these requirements.

The Model 335 protects the operator and surrounding area from electric shock or burn, mechanical hazards, excessive temperature, and spread of fire from the instrument. Environmental conditions outside of the conditions below may pose a hazard to the operator and surrounding area.

- Indoor use
- Altitude to 2000 m
- Temperature for safe operation: 5 °C to 40 °C
- Maximum relative humidity: 80% for temperature up to 31 °C decreasing linearly to 50% at 40 °C
- Power supply voltage fluctuations not to exceed  $\pm 10\%$  of the nominal voltage
- Overvoltage category II
- Pollution degree 2
- IPX0: not protected against harmful ingress of water

### *Ground the Instrument*

To minimize shock hazard, the instrument is equipped with a 3-conductor AC power cable. Plug the power cable into an approved 3-contact electrical outlet or use a 3-contact adapter with the grounding wire (green) firmly connected to an electrical ground (safety ground) at the power outlet. The power jack and mating plug of the power cable meet Underwriters Laboratories (UL) and International Electrotechnical Commission (IEC) safety standards.

### *Ventilation*

The instrument has ventilation holes in its side covers. Do not block these holes when the instrument is operating.

### *Do Not Operate in an Explosive Atmosphere*

Do not operate the instrument in the presence of flammable gases or fumes. Operation of any electrical instrument in such an environment constitutes a definite safety hazard.

### *Keep Away from Live Circuits*

Operating personnel must not remove instrument covers. Refer component replacement and internal adjustments to qualified maintenance personnel. Do not replace components with power cable connected. To avoid injuries, always disconnect power and discharge circuits before touching them.

### *Do Not Substitute Parts or Modify Instrument*

Do not install substitute parts or perform any unauthorized modification to the instrument. Return the instrument to an authorized Lake Shore Cryotronics, Inc. representative for service and repair to ensure that safety features are maintained.

### *Cleaning*

Do not submerge instrument. Clean only with a damp cloth and mild detergent. Exterior only.













	Direct current (power line)		Equipment protected throughout by double insulation or reinforces insulation (equivalent to Class II of IEC 536—see Annex H)
	Alternating current (power line)		
	Alternating or direct current (power line)		CAUTION: High voltages; danger of electric shock; background color: yellow; symbol and outline: black
	Three-phase alternating current (power line)		
	Earth (ground) terminal		
	Protective conductor terminal		CAUTION or WARNING: See instrument documentation; background color: yellow; symbol and outline: black
	Frame or chassis terminal		
	On (supply)		
	Off (supply)		

FIGURE 1-4 *Safety symbols*

## 6.3 USB Interface

The Model 335 USB interface provides a convenient way to connect to most modern computers, as a USB interface is provided on nearly all new PCs as of the writing of this manual. The USB interface is implemented as a virtual serial com port connection. This implementation provides a simple migration path for modifying existing RS-232 based remote interface software. It also provides a simpler means of communicating than a standard USB implementation.

### 6.3.1 Physical Connection

The Model 335 has a B-type USB connector on the rear panel. This is the standard connector used on USB peripheral devices, and it allows the common USB A-type to B-type cable to be used to connect the Model 335 to a host PC. The pin assignments for A-type and B-type connectors are shown in section 8.10. The maximum length of a USB cable, as defined by the USB 2.0 standard, is 5 m (16.4 ft). This length can be extended using USB hubs every 5 m (16.4 ft) up to 5 times, for a maximum total length of 30 m (98.4 ft).

### 6.3.2 Hardware Support

The USB interface emulates an RS-232 serial port at a fixed 57,600 baud rate, but with the physical connections of a USB. This programming interface requires a certain configuration to communicate properly with the Model 335. The proper configuration parameters are listed in TABLE 6-5.

<b>Baud rate</b>	57,600
<b>Data bits</b>	7
<b>Start bits</b>	1
<b>Stop bits</b>	1
<b>Parity</b>	Odd
<b>Flow control</b>	None
<b>Handshaking</b>	None

TABLE 6-5 *Host com port configuration*

The USB hardware connection uses the full speed (12,000,000 bits/s) profile of the USB 2.0 standard; however, since the interface uses a virtual serial com port at a fixed data rate, the data throughput is still limited to a baud rate of 57,600 bits/s.

### 6.3.3 Installing the USB Driver

The Model 335 USB driver has been made available through Windows® Update. This is the recommended method for installing the driver, as it will ensure that you always have the latest version of the driver installed. If you are unable to install the driver from Windows® Update, refer to section 6.3.3.3 to install the driver from the web.

These procedures assume that you are logged into a user account that has administrator privileges.

#### 6.3.3.1 Installing the Driver From Windows® Update in Windows® 7 and Vista®

1. Connect the USB cable from the Model 335 to the computer.
2. Turn on the Model 335.
3. When the Found New Hardware wizard appears, select **Locate and install driver software (recommended)**.
4. If User Account Control (UAC) is enabled, a UAC dialog box may appear asking if you want to continue. Click **Continue**.
5. The Found New Hardware wizard should automatically connect to Windows® Update and install the drivers.



If the Found New Hardware wizard is unable to connect to Windows® Update or find the drivers, a message to “Insert the disc that came with your Lake Shore Model 335” will be displayed. Click Cancel and refer to section 6.3.3.3 to install the driver from the web.

6. When the Found New Hardware wizard finishes installing the driver, a confirmation message stating “the software for this device has been successfully installed” will appear. Click Close to complete the installation.

#### 6.3.3.2 Installing the Driver From Windows® Update in Windows® XP

1. Connect the USB cable from the Model 335 to the computer.
2. Turn on the Model 335.
3. When the Found New Hardware wizard appears, select **Yes, this time only** and click **Next**.
4. Select **Install the software automatically (Recommended)** and click **Next**.
5. The Found New Hardware wizard should automatically connect to Windows® Update and install the drivers.



If the Found New Hardware wizard is unable to connect to Windows® Update or find the drivers, a message saying **Cannot Install this Hardware** will be displayed. Click the **Cancel** button and refer to section 6.3.3.3 to install the driver from the web.

6. When the Found New Hardware wizard finishes installing the driver a confirmation message stating “the wizard has finished installing the software for Lake Shore Model 335 Temperature Controller” will appear. Click **Finish** to complete the installation.

#### 6.3.3.3 Installing the Driver From the Web

The Model 335 USB driver is available on the Lake Shore website. To install the driver it must be downloaded from the website and extracted. Use the procedure in section 6.3.3.1 through section 6.3.4 to download, extract, and install the driver using Windows® 7, Vista® and XP.

##### 6.3.3.3.1 Download the driver:

1. Locate the Model 335 USB driver on the downloads page on the Lake Shore website.
2. Right-click on the USB driver download link, and select **save target/link as**.
3. Save the driver to a convenient place, and take note as to where the driver was downloaded.

##### 6.3.3.3.2 Extract the driver:

The downloaded driver is in a ZIP compressed archive. The driver must be extracted from this file. Windows® provides built-in support for ZIP archives. If this support is disabled, a third-party application, such as WinZip™ or 7-Zip, must be used.

For Windows® 7 and Vista®:

1. Right click on the file and click **extract all**.
2. An Extract Compressed (Zipped) Folders dialog box will appear. It is recommended the default folder is not changed. Take note of this folder location.
3. Click to clear the **Show extracted files when complete** checkbox, and click **Extract**.

For Windows® XP:

1. Right-click on the file and click **extract all**.
2. The Extraction wizard will appear. Click **Next**.

3. It is recommended the default folder is not changed. Take note of this folder location and click **Next**.
4. An "Extraction complete" message will be displayed. Click to clear the **Show extracted files** checkbox, and click **Finish**.

#### 6.3.3.3 Manually install the driver

Manually installing drivers differ between versions of Windows®. The following sections describe how to manually install the driver using Windows® 7, Vista® and XP. To install the driver you must be logged into a user account that has administrator privileges.

For Windows® 7 and Vista®:

1. Connect the USB cable from the Model 335 to the computer.
2. Turn on the Model 335.
3. If the Found New Hardware wizard appears, click **Ask me again later**.
4. Open Device Manager. Use this procedure to open Device Manager.
  - a. Click the Windows® **Start** button and type Device Manager in the **Start Search** box.
  - b. Click on the Device Manager link in the Search Results Under Programs dialog box.
  - c. If User Account Control is enabled click **Continue** on the User Account Control prompt.
5. Click **View** and ensure the **Devices by Type** check box is selected.
6. In the main window of Device Manager, locate **Other Devices** in the list of device types. In many instances this will be between Network adapters and Ports (COM & LPT). If the **Other Devices** item is not already expanded, click the + icon. Lake Shore Model 335 should appear indented underneath **Other Devices**. If it is not displayed as Lake Shore Model 335, it might be displayed as USB Device. If neither are displayed, click **Action** and then **Scan for hardware changes**, which may open the Found New Hardware wizard automatically. If the Found New Hardware wizard opens, click **Cancel**.
7. Right-click on Lake Shore Model 335 and click **Update Driver Software**.
8. Click **Browse my computer for driver software**.
9. Click **Browse** and select the location of the extracted driver.
10. Ensure the **Include subfolders** check box is selected and click **Next**.
11. When the driver finishes installing a confirmation message stating "Windows has successfully updated your driver software" should appear. Click **Close** to complete the installation.

For Windows® XP:

1. Connect the USB cable from the Model 335 to the computer.
2. Turn on the Model 335.
3. The Found New Hardware wizard should appear. If the Found New Hardware wizard does not appear, the following procedure can be used to open the Hardware Update wizard which can be used instead:
  - a. Open Device Manager. Use this procedure to open the Device Manager:
    - Right-click on **My Computer** and then click **Properties**. This will open the System Properties dialog.
    - Click the **Hardware** tab and then click **Device Manager**.
  - b. Click **View** and ensure the **Devices by Type** check box is selected.

- c. In the main window of Device Manager, locate the **Ports (COM & LPT)** device type. In many instances this will be between the Network adapters and Processors items. If the **Ports (COM & LPT)** item is not already expanded, click the + icon. Lake Shore Model 335 should appear indented underneath **Ports (COM & LPT)**. If it is not displayed as Lake Shore Model 335, it might be displayed as USB Device. If neither are displayed, click **Action** and then select **Scan for hardware changes**, which may open the Found New Hardware wizard automatically. If the Found New Hardware wizard opens, continue to step 4.
- d. Right-click on Lake Shore Model 335 and click **Update Driver**.
4. Select **No, not at this time** and click **Next**.
5. Select **Search for the best driver in these locations**, click to clear the **Search removable media** check box, and click the **Include this location in the search** check box.
6. Click **Browse** and open the location of the extracted driver.
7. Click **Next**.
8. When the driver finishes installing a confirmation message stating “The wizard has finished installing the software for Lake Shore Model 335 Temperature Controller” should appear. Click **Finish** to complete the installation.

## 6.3.4 Communication

Communicating via the USB interface is done using message strings. The message strings should be carefully formulated by the user program according to some simple rules to establish effective message flow control.

### 6.3.4.1 Character Format

A character is the smallest piece of information that can be transmitted by the interface. Each character is ten bits long and contains data bits, bits for character timing, and an error detection bit. The instrument uses 7 bits for data in the American Standard Code for Information Interchange (ASCII) format. One start bit and one stop bit are necessary to synchronize consecutive characters. Parity is a method of error detection. One parity bit configured for odd parity is included in each character.

ASCII letter and number characters are used most often as character data. Punctuation characters are used as delimiters to separate different commands or pieces of data. A special ASCII character, line feed (LF OAH), is used to indicate the end of a message string. This is called the message terminator. The Model 335 will accept either the line feed character alone, or a carriage return (CR ODH) followed by a line feed as the message terminator. The instrument query response terminator will include both carriage return and line feed.

### 6.3.4.2 Message Strings

A message string is a group of characters assembled to perform an interface function. There are three types of message strings: commands, queries and responses. The computer issues command and query strings through user programs, the instrument issues responses. Two or more command or query strings can be chained together in one communication, but they must be separated by a semi-colon (;). The total communication string must not exceed 255 characters in length.

A command string is issued by the computer and instructs the instrument to perform a function or change a parameter setting. The format is:

**<command mnemonic><space><parameter data><terminators>.**

Command mnemonics and parameter data necessary for each one is described in section 6.4. Terminators must be sent with every message string.



A query string is issued by the computer and instructs the instrument to send a response. The query format is:

**<query mnemonic><?><space><parameter data><terminators>.**

Query mnemonics are often the same as commands with the addition of a question mark. Parameter data is often unnecessary when sending queries. Query mnemonics and parameter data if necessary is described in section 6.4. Terminators must be sent with every message string. The computer should expect a response very soon after a query is sent.

A response string is the instrument's response or answer to a query string. The response can be a reading value, status report or the present value of a parameter. Response data formats are listed along with the associated queries in section 6.4. The response is sent as soon as possible after the instrument receives the query.

### 6.3.5 Message Flow Control

It is important to remember that the user program is in charge of the USB communication at all times. The instrument cannot initiate communication, determine which device should be transmitting at a given time, or guarantee timing between messages. All of this is the responsibility of the user program.

When issuing commands the user program alone should:

- Properly format and transmit the command including the terminator as one string
- Guarantee that no other communication is started for 50 ms after the last character is transmitted
- Not initiate communication more than 20 times/s

When issuing queries or queries and commands together, the user program should:

- Properly format and transmit the query including the terminator as one string
- Prepare to receive a response immediately
- Receive the entire response from the instrument including the terminator
- Guarantee that no other communication is started during the response or for 50 ms after it completes
- Not initiate communication more than 20 times/s

Failure to follow these simple rules will result in inability to establish communication with the instrument or intermittent failures in communication.

6.4 Command Summary

This section provides a listing of the interface commands. A summary of all the commands is provided in TABLE 6-6. All the commands are detailed in section 6.4.1, and are presented in alphabetical order.




FIGURE 6-5 Sample command format




FIGURE 6-6 Sample query format

Command	Function	Page	Command	Function	Page
*CLS	Clear Interface Cmd	107	INNAME	Sensor Input Name Cmd	117
*ESE	Event Status Enable Register Cmd	107	INNAME?	Sensor Input Name Query	117
*ESE?	Event Status Enable Register Query	107	INTYPE	Input Type Parameter Cmd	118
*ESR?	Standard Event Status Register Query	108	INTYPE?	Input Type Parameter Query	118
*IDN?	Identification Query	108	KRDG?	Kelvin Reading Query	119
*OPC	Operation Complete Cmd	108	LEDS	Front Panel LEDS Cmd	119
*OPC?	Operation Complete Query	108	LEDS?	Front Panel LEDS Query	119
*RST	Reset Instrument Cmd	108	LOCK	Front Panel Keyboard Lock Cmd	119
*SRE	Service Request Enable Register Cmd	109	LOCK?	Front Panel Keyboard Lock Query	119
*SRE?	Service Request Enable Register Query	109	MDAT?	Minimum/Maximum Data Query	119
*STB?	Status Byte Query	109	MNMXRST	Minimum and Maximum Function Reset Cmd	119
*TST?	Self-Test Query	109	MODE	Remote Interface Mode Cmd	120
*WAI	Wait-to-Continue Cmd	109	MODE?	Remote Interface Mode Query	120
ALARM	Input Alarm Parameter Cmd	110	MOUT	Manual Output Cmd	120
ALARM?	Input Alarm Parameter Query	110	MOUT?	Output Manual Heater Power (MHP) Output Query	120
ALARMST?	Input Alarm Status Query	110	OPST?	Operational Status Query	120
ALMRST	Reset Alarm Status Cmd	110	OPSTE	Operation Status Enable Cmd	120
ANALOG	Monitor Out Parameter Cmd	111	OPSTE?	Operational Status Enable Query	120
ANALOG?	Monitor Out Parameter Query	111	OPSTR?	Operational Status Register Query	121
ATUNE	Autotune Cmd	111	OUTMODE	Output Mode Command	121
BRIGT	Display Brightness Cmd	111	OUTMODE?	Output Mode Query	121
BRIGT?	Display Brightness Query	112	POLARITY	Output Voltage Polarity Command	122
CRDG?	Celsius Reading Query	112	POLARITY?	Output Voltage Polarity Query	122
CRVDEL	Curve Delete Cmd	112	PID	Control Loop PID Values Cmd	121
CRVHDR	Curve Header Cmd	112	PID?	Control Loop PID Values Query	121
CRVHDR?	Curve Header Query	112	RAMP	Control Setpoint Ramp Parameter Cmd	122
CRVPT	Curve Data Point Cmd	113	RAMP?	Control Setpoint Ramp Parameter Query	122
CRVPT?	Curve Data Point Query	113	RAMPST?	Control Setpoint Ramp Status Query	122
DFLT	Factory Defaults Cmd	113	RANGE	Heater Range Cmd	123
DIOCUR	Diode Excitation Current Parameter Cmd	113	RANGE?	Heater Range Query	123
DIOCUR?	Diode Excitation Current Parameter Query	113	RDGST?	Input Reading Status Query	123
DISPFLD	Custom ModeDisplay Field Cmd	114	RELAY	Relay Control Parameter Cmd	123
DISPFLD?	Custom Mode Display Field Query	114	RELAY?	Relay Control Parameter Query	123
DISPLAY	Display Setup Cmd	114	RELAYST?	Relay Status Query	124
DISPLAY?	Display Setup Query	114	SCAL	Generate SoftCal Curve Cmd	124
EMUL	Model 331/332 Emulation Mode Cmd	115	SETP	Control Setpoint Cmd	124
EMUL?	Model 331/332 Emulation Mode Query	115	SETP?	Control Setpoint Query	124
FILTER	Input Filter Parameter Cmd	115	SRDG?	Sensor Units Input Reading Query	125
FILTER?	Input Filter Parameter Query	115	TEMP?	Thermocouple Junction Temperature Query	125
HTR?	Heater Output Query	115	TLIMIT	Temperature Limit Cmd	125
HTRSET	Heater Setup Cmd	116	TLIMIT?	Temperature Limit Query	125
HTRSET?	Heater Setup Query	116	TUNEST?	Control Tuning Status Query	125
HTRST?	Heater Status Query	116	WARMUP	Warmup Supply Parameter Cmd	126
IEEE	IEEE-488 Parameter Cmd	116	WARMUP?	Warmup Supply Parameter Query	126
IEEE?	IEEE-488 Interface Parameter Query	117	ZONE	Control Loop Zone Table Parameter Cmd	126
INCRV	Input Curve Number Cmd	117	ZONE?	Output Zone Table Parameter Query	126
INCRV?	Input Curve Number Query	117			

TABLE 6-6 Command summary

## 6.4.1 Interface Commands

This section lists the interface commands in alphabetical order.

<b>*</b>	Begins common interface command
<b>?</b>	Required to identify queries
<b>s[n]</b>	String of alphanumeric characters with length “n.” Send these strings using surrounding quotes. Quotes enable characters such as commas and spaces to be used without the instrument interpreting them as delimiters.
<b>nn...</b>	String of number characters that may include a decimal point.
<b>dd</b>	Dotted decimal format, common with IP addresses. Always contains 4 dot separated 3-digit decimal numbers, such as 192.168.000.012.
<b>[term]</b>	Terminator characters
<b>&lt;...&gt;</b>	Indicates a parameter field, many are command specific.
<b>&lt;state&gt;</b>	Parameter field with only On/Off or Enable/Disable states.
<b>&lt;value&gt;</b>	Floating point values have varying resolution depending on the type of command or query issued.

TABLE 6-7 Interface commands key

### \*CLS

Input

Remarks

### Clear Interface Command

**\*CLS [term]**

Clears the bits in the Status Byte Register, Standard Event Status Register, and Operation Event Register, and terminates all pending operations. Clears the interface, but not the controller. The related controller command is **\*RST**.

### \*ESE

Input

Format

Remarks

### Event Status Enable Register Command

**\*ESE <bit weighting> [term]**

nnn

Each bit is assigned a bit weighting and represents the enable/disable mask of the corresponding event flag bit in the Standard Event Status Register. To enable an event flag bit, send the command **\*ESE** with the sum of the bit weighting for each desired bit. Refer to section 6.2.5 for a list of event flags.

Example

To enable event flags 0, 4, and 7, send the command **\*ESE 145[term]**. 145 is the sum of the bit weighting for each bit.

Bit	Bit Weighting	Event Name
0	1	OPC
2	4	QXE
4	16	EXE
5	32	CME
7	128	PON
Total:	181	

### \*ESE?

Input

Returned

Format

### Event Status Enable Register Query

**\*ESE? [term]**

<bit weighting>[term]

nnn (refer to section 6.2.5 for a list of event flags)

<hr/>	
<b>*ESR?</b>	<b>Standard Event Status Register Query</b>
Input	<b>*ESR?</b> [ <i>term</i> ]
Returned	<bit weighting>
Format	nnn
Remarks	The integer returned represents the sum of the bit weighting of the event flag bits in the Standard Event Status Register. Refer to section 6.2.5 for a list of event flags.
<hr/>	
<b>*IDN?</b>	<b>Identification Query</b>
Input	<b>*IDN?</b> [ <i>term</i> ]
Returned	<manufacturer>,<model>,<instrument serial>/<option serial>,<firmware version>[ <i>term</i> ]
Format	s[4],s[8],s[7]/s[7],n.n <manufacturer>      Manufacturer ID <model>            Instrument model number <instrument serial>   Instrument serial number <option card serial>   Option card serial number <firmware version>   Instrument firmware version
Example	LSCI,MODEL335,1234567/1234567,1.0
<hr/>	
<b>*OPC</b>	<b>Operation Complete Command</b>
Input	<b>*OPC</b> [ <i>term</i> ]
Remarks	Generates an Operation Complete event in the Event Status Register upon completion of all pending selected device operations. Send it as the last command in a command string.
<hr/>	
<b>*OPC?</b>	<b>Operation Complete Query</b>
	<b>*OPC?</b> [ <i>term</i> ]
Returned	1[ <i>term</i> ]
Remarks	Places a 1 in the controller output queue upon completion of all pending selected device operations. Send as the last command in a command string. Not the same as *OPC.
<hr/>	
<b>*RST</b>	<b>Reset Instrument Command</b>
Input	<b>*RST</b> [ <i>term</i> ]
Remarks	Sets controller parameters to power-up settings.

---

**\*SRE Service Request Enable Register Command**

Input	<b>*SRE</b> <bit weighting>[term]
Format	nnn
Remarks	Each bit has a bit weighting and represents the enable/disable mask of the corresponding status flag bit in the Status Byte Register. To enable a status flag bit, send the command <b>*SRE</b> with the sum of the bit weighting for each desired bit. Refer to section 6.2.6 for a list of status flags.
Example	To enable status flags 4, 5, and 7, send the command <b>*SRE 208</b> [term]. 208 is the sum of the bit weighting for each bit.

Bit	Bit Weighting	Event Name
4	16	MAV
5	64	ESB
7	128	OSB
Total:	208	

---

**\*SRE? Service Request Enable Register Query**

Input	<b>*SRE?</b> [term]
Returned	<bit weighting>[term]
Format	nnn (refer to section 6.2.6 for a list of status flags)

---

**\*STB? Status Byte Query**

Input	<b>*STB?</b> [term]
Returned	<bit weighting>[term]
Format	nnn
Remarks	Acts like a serial poll, but does not reset the register to all zeros. The integer returned represents the sum of the bit weighting of the status flag bits that are set in the Status Byte Register. Refer to section 6.2.6 for a list of status flags.

---

**\*TST? Self-Test Query**

Input	<b>*TST?</b> [term]
Returned	<status>[term]
Format	n <status>      0 = no errors found, 1 = errors found
Remarks	The Model 335 reports status based on the test done at power up.

---

**\*WAI Wait-to-Continue Command**

Input	<b>*WAI</b> [term]
Remarks	Causes the IEEE-488 interface to hold off until all pending operations have been completed. This is the same function as the <b>*OPC</b> command, except that it does not set the Operation Complete event bit in the Event Status Register.

---

**ALARM**      **Input Alarm Parameter Command**

Input	<b>ALARM</b> <input>,<off/on>,<high value>,<low value>,<deadband>,<latch enable>,<audible>,<visible> [term]
Format	a,n,±nnnnnn,±nnnnnn,±nnnnnn,n,n,n <input>            Specifies which input to configure: A or B. <off/on>           Determines whether the instrument checks the alarm for this input, where 0 = off and 1 = on. <high setpoint>   Sets the value the source is checked against to activate the high alarm. <low setpoint>    Sets the value the source is checked against to activate low alarm. <deadband>        Sets the value that the source must change outside of an alarm condition to deactivate an unlatched alarm. <latch enable>     Specifies a latched alarm (remains active after alarm condition correction) where 0 = off (no latch) and 1 = on. <audible>          Specifies if the internal speaker will beep when an alarm condition occurs. Valid entries: 0 = off, 1 = on. <visible>          Specifies if the Alarm LED on the instrument front panel will blink when an alarm condition occurs. Valid entries: 0 = off, 1 = on
Remarks	Configures the alarm parameters for an input.
Example	<b>ALARM A,0[term]</b> turns off alarm checking for Input A. <b>ALARM B,1,270.0,0,0,1,1,1[term]</b> turns on alarm checking for input B, activates high alarm if kelvin reading is over 270, and latches the alarm when kelvin reading falls below 270. Alarm condition will cause instrument to beep and the front panel Alarm LED to blink.

---

**ALARM?**      **Input Alarm Parameter Query**

Input	<b>ALARM?</b> <input> [term]
Format	a <input>            A or B
Returned	<off/on>,<high value>,<low value>,<deadband>,<latch enable>,<audible>,<visible> [term]
Format	n,±nnnnnn,±nnnnnn,±nnnnnn,n,n,n (refer to command for description)

---

**ALARMST?**    **Input Alarm Status Query**

Input	<b>ALARMST?</b> <input> [term]
Format	a <input>            A - B
Returned	<high state>,<low state>[term]
Format	n,n <high state> 0 = Off, 1 = On <low state> 0 = Off, 1 = On

---

**ALMRST**      **Reset Alarm Status Command**

Input	<b>ALMRST</b> [term]
Remarks	Clears both the high and low status of all alarms, including latching alarms.

---

**ANALOG Monitor Out Parameter Command****Input** **ANALOG** <output>,<input>,<units>,<high value>,<low value>,<polarity>[term]**Format** n,n,n,±nnnnn,±nnnnn,n

<output> Output 2 is the only valid entry and must be included (included for compatibility with other Lake Shore temperature instruments). 2 = Output 2

<input> Specifies which input to monitor. 0 = none, 1 = Input A, 2 = Input B

<units> Specifies the units on which to base the output voltage: 1 = kelvin, 2 = Celsius, 3 = sensor units

<high value> If output mode is Monitor Out, this parameter represents the data at which the Monitor Out reaches +100% output. Entered in the units designated by the <units> parameter. Refer to OUTMODE command.

<low value> If output mode is Monitor Out, this parameter represents the data at which the analog output reaches -100% output if bipolar, or 0% output if positive only. Entered in the units designated by the <units> parameter.

<polarity> Specifies output voltage is 0 = unipolar (positive output only) or 1 = bipolar (positive or negative output)

**Example** **ANALOG 2,1,1,100.0,0.0,0[term]** sets Output 2 Voltage mode to monitor Input A kelvin reading with 100.0 K at +100% output (+10.0 V) and 0.0 K at 0% output (0.0 V).**Remarks** Use the OUTMODE command to set the output mode to Monitor Out. The <input> parameter in the ANALOG command is the same as the <input> parameter in the OUTMODE command. It is included in the ANALOG command for backward compatibility with previous Lake Shore temperature monitors and controllers. The ANALOG command name is also named as such for backward compatibility.

---

**ANALOG? Monitor Out Parameter Query****Input** **ANALOG?** <output>[term]**Format** n**Returned** <input>,<units>,<high value>,<low value>,<polarity>[term]**Format** n,n,±nnnnn,±nnnnn,n (refer to command for definition)

---

**ATUNE Autotune Command****Input** **ATUNE** <output>,<mode>,[term]**Format** n,n

&lt;output&gt; Specifies the output associated with the loop to be Autotuned: 1 or 2.

&lt;mode&gt; Specifies the Autotune mode. Valid entries: 0 = P Only, 1 = P and I, 2 = P, I, and D.

**Example** **ATUNE 2,1[term]** initiates Autotuning of control loop associated with output 2, in P and I mode.**Remarks** If initial conditions required to Autotune the specified loop are not met, an Autotune initialization error will occur and the Autotune process will not be performed. The TUNEST? query can be used to check if an Autotune error occurred.

---

**BRIGT Display Brightness Command****Input** **BRIGT** <brightness value>[term]**Format** n

&lt;brightness value&gt; 0–3

**Remarks** Sets the display brightness for the front panel display 0=25%, 1=50%, 2=75%, 3=100%.



<b>BRIGT?</b> Input Returned Format	<b>Display Brightness Query</b> <b>BRIGT?</b> [ term] <brightness value>[term] n (refer to command for description)
<b>CRDG?</b> Input Format Returned Format Remarks	<b>Celsius Reading Query</b> <b>CRDG?</b> <input> [ term] a <input> A or B <temp value>[term] ±nnnnnn Also see the RDGST? query.
<b>CRVDEL</b> Input Format Example	<b>Curve Delete Command</b> <b>CRVDEL</b> <curve> [ term] nn <curve> Specifies a user curve to delete. Valid entries: 21–59. <b>CRVDEL 21[term]</b> deletes User Curve 21.
<b>CRVHDR</b> Input Format Remarks Example	<b>Curve Header Command</b> <b>CRVHDR</b> <curve>, <name>, <SN>, <format>, <limit value>, <coefficient> [ term] nn,s[15],s[10],n,+nnn.nnn,n <curve> Specifies which curve to configure. Valid entries: 21–59. <name> Specifies curve name. Limited to 15 characters. <SN> Specifies the curve serial number. Limited to ten characters. <format> Specifies the curve data format. Valid entries: 1 = mV/K, 2 = V/K, 3 = Ω/K, 4 = log Ω/K. <limit value> Specifies the curve temperature limit in kelvin. <coefficient> Specifies the curves temperature coefficient. Valid entries: 1 = negative, 2 = positive. Configures the user curve header. The coefficient parameter will be calculated automatically based on the first two curve datapoints. It is included as a parameter for compatability with the CRVHDR? query. All parameter fields must be filled with a value. If the sensor does not have a serial number, enter NONE in this parameter field. <b>CRVHDR 21,DT-470,00011134,2,325.0,1[term]</b> configures User Curve 21 with a name of DT-470, serial number of 00011134, data format of volts versus kelvin, upper temperature limit of 325 K, and negative coefficient.
<b>CRVHDR?</b> Input Format Returned Format	<b>Curve Header Query</b> <b>CRVHDR?</b> <curve> [ term] nn <curve> Valid entries: 1–59. <name>,<SN>,<format>,<limit value>,<coefficient>[term] s[15],s[10],n,+nnn.nnn,n (refer to command for description)

---

**CRVPT Curve Data Point Command**

Input	<b>CRVPT</b> <curve>,<index>,<units value>,<temp value>[term]
Format	nn,nnn,±nnnnnn,±nnnnnn
	<curve> Specifies which curve to configure. Valid entries: 21–59.
	<index> Specifies the points index in the curve. Valid entries: 1–200.
	<units value> Specifies sensor units for this point to six digits.
	<temp value> Specifies the corresponding temperature in kelvin for this point to six digits.
Remarks	Configures a user curve data point.
Example	<b>CRVPT 21,2,0.10191,470.000,N[term]</b> sets User Curve 21 second data point to 0.10191 sensor units and 470.000 K.

---

**CRVPT? Curve Data Point Query**

Input	<b>CRVPT?</b> <curve>,<index>[term]
Format	nn,nnn
	<curve> Specifies which curve to query: 1–59.
	<index> Specifies the points index in the curve: 1–200.
Returned	<units value>,<temp value>[term]
Format	±nnnnnn,±nnnnnn (refer to command for description)
Remarks	Returns a standard or user curve data point.

---

**DFLT Factory Defaults Command**

Input	<b>DFLT 99</b> [term]
Remarks	Sets all configuration values to factory defaults and resets the instrument. The “99” is included to prevent accidentally setting the unit to defaults.

---

**DIOCUR Diode Excitation Current Parameter Command**

Input	<b>DIOCUR</b> <input>,<excitation>[term]
Format	a,n
	<input> Specifies which input to configure: A or B.
	<excitation> Specifies the Diode excitation current: 0 = 10 $\mu$ A, 1 = 1 mA.
Remarks	The 10 $\mu$ A excitation current is the only calibrated excitation current, and is used in almost all applications. Therefore the Model 335 will default the 10 $\mu$ A current setting any time the input sensor type is changed in order to prevent an accidental change. If using a current that is not 10 $\mu$ A, the input sensor type must first be configured to Diode (INTYPE command). If the sensor type is not set to Diode when the DIOCUR command is sent, the command will be ignored.

---

**DIOCUR? Diode Excitation Current Parameter Query**

Input	<b>DIOCUR?</b> <input>[term]
Format	a
	<input> A or B
Returned	<excitation>[term]
Format	n (refer to command for description)

<b>DISPFLD</b>	<b>Custom Mode Display Field Command</b>
Input Format	<b>DISPFLD</b> <field>,<source>,<units>[term] n,n,n <field> Specifies field (display location) to configure: 1–4. <source> Specifies item to display in the field: 0 = None, 1 = Input A, 2 = Input B, 3 = Setpoint 1, 4 = Setpoint 2, 5=Output 1, 6=Output 2 <units> Valid entries: 1 = kelvin, 2 = Celsius, 3 = sensor units, 4 = minimum data, 5 = maximum data, 6 = sensor name.
<b>Example</b>	<b>DISPFLD 2,1,1[term]</b> displays kelvin reading for Input A in display field 2 when display mode is set to Custom.
<b>Remarks</b>	Since each display field is only ten characters, only the first nine characters of the sensor name will be displayed when the units for a field are set to display the sensor name. However, if two adjacent fields are assigned to the same sensor name, then the entire twenty character line can be used, allowing all fifteen sensor name characters to be displayed.
This command only applies to the readings displayed in the Custom display mode. All other display modes have predefined readings in predefined locations, and will use the Preferred Units parameter to determine which units to display for each sensor input. Refer to section 4.3 for details on display setup	
<b>DISPFLD?</b>	<b>Custom Mode Display Field Query</b>
Input Format	<b>DISPFLD?</b> <field>[term] n <field> Specifies field (display location) to query: 1–4.
Returned Format	<input>,<units>[term] n,n (refer to command for description)
<b>DISPLAY</b>	<b>Display Setup Command</b>
Input Format	<b>DISPLAY</b> <mode>[term] n <mode> Specifies display mode: 0 = Input A, 1 = Input A Max/Min, 2 = Two Input, Loop A, 3=Input B, 4=Input B Max/Min, 5=Two Input, Loop B, 6=Custom, 7=Two Loop
<b>Remarks</b>	When the input display mode is set to Custom, use the DISPFLD command to configure the display.
<b>DISPLAY?</b>	<b>Display Setup Query</b>
Input Returned Format	<b>DISPLAY?</b> [term] <mode>[term] n (refer to command for description)

<b>EMUL</b> Input Format	<b>Model 331/332 Interface Emulation Mode Command</b> <b>EMUL</b> <emulation mode>,<PID scaling mode>[term] nn <emulation mode>      0=None (335), 1=331, 2=332 <PID scaling mode>      PID control scaling: 0=335 (Temperature), 1=331/332 (Sensor)
Remarks	The 331 and 332 emulation modes provide a means of using the Model 335 in place of a Model 331 or 332 in a software controlled system without updating the software. The emulation mode setting only affects remote operation; front panel operation of the Model 335 is not changed. Please reference the Model 331 or 332 user's manual for information on the interface commands for those models. For more information on using the emulation modes, see section 5.11.
<b>EMUL?</b> Input Returned Format	<b>Model 331/332 Interface Emulation Mode Query</b> <b>EMUL?</b> [term] <emulation mode>,<PID scaling mode>[term] n,n (refer to command for description)
<b>FILTER</b> Input Format	<b>Input Filter Parameter Command</b> <b>FILTER</b> <input>,<off/on>,<points>,<window>[term] a,n,nn,nn <input>      Specifies input to configure: A or B. <off/on>      Specifies whether the filter function is 0 = Off or 1 = On. <points>      Specifies how many data points the filtering function uses. Valid range = 2–64. <window>      Specifies what percent of full scale reading limits the filtering function. Reading changes greater than this percentage reset the filter. Valid range = 1 to 10%.
Example	<b>FILTER B,1,10,2[term]</b> filter input B data through ten readings with 2% of full scale window.
<b>FILTER?</b> Input Format  Returned Format	<b>Input Filter Parameter Query</b> <b>FILTER?</b> <input>[term] a <input>      Specifies input to query: A or B. <off/on>,<points>,<window>[term] n,nn,nn (refer to command for description)
<b>HTR?</b> Input Format  Returned Format	<b>Heater Output Query</b> <b>HTR?</b> <output>[term] n <output>      Heater output to query: 1 = Output 1, 2 = Output 2 <heater value>[term] ±nnn.n <heater value>      Percent of full scale current for Output 1 and Output 2 in Current mode, or percent of full scale voltage for Output 2 in Voltage mode.

<hr/>	
<b>HTRSET</b>	<b>Heater Setup Command</b>
Input	<b>HTRSET</b> <output>,<type>,<heater resistance>,<max current>,<max user current>,<current/power>[term]
Format	n,n,n,n,+n.nnn,n
	<div> <div>&lt;output&gt;</div> <div>Specifies which heater output to configure: 1 or 2.</div> </div> <div> <div>&lt;type&gt;</div> <div>Output type (Output 2 only): 0=Current, 1=Voltage</div> </div> <div> <div>&lt;htr resistance&gt;</div> <div>Heater Resistance Setting: 1 = 25 <math>\Omega</math>, 2 = 50 <math>\Omega</math>.</div> </div> <div> <div>&lt;max current&gt;</div> <div>Specifies the maximum heater output current: 0 = User Specified, 1 = 0.707 A, 2 = 1 A, 3 = 1.141 A, 4 = 1.732 A</div> </div> <div> <div>&lt;max user current&gt;</div> <div>Specifies the maximum heater output current if max current is set to User Specified.</div> </div> <div> <div>&lt;current/power&gt;</div> <div>Specifies whether the heater output displays in current or power (current mode only). Valid entries: 1 = current, 2 = power.</div> </div>
Example	<b>HTRSET 2,0,1,2,0,1[term]</b> Heater Output 2 uses the current source output, the 25 $\Omega$ heater setting, has a maximum current of 1 A, the maximum user current is set to 0 A because it is not going to be used since a discrete value has been chosen, and the heater output will be displayed in units of current.
Remarks	Max current will be limited to 1 A on Output 1 if the heater resistance is set to 50 $\Omega$ . If the heater resistance is set to 25 $\Omega$ on Output 1, then the max current will be limited to 1.414 A if Output 2 is in Current mode, or 1.732 A if Output 2 is in Voltage mode. See section 4.5.1.3 for more information on the max current setting.
<hr/>	
<b>HTRSET?</b>	<b>Heater Setup Query</b>
Input	<b>HTRSET?</b> <output>[term]
Format	n
	<div> <div>&lt;output&gt;</div> <div>Specifies which heater output to query: 1 or 2.</div> </div>
Returned	<type>,<htr resistance>,<max current>,<max user current>,<current/power>[term]
Format	n,n,n,+n.nnn,n
<hr/>	
<b>HTRST?</b>	<b>Heater Status Query</b>
Input	<b>HTRST?</b> <output>[term]
Format	n
	<div> <div>&lt;output&gt;</div> <div>Specifies which heater output to query: 1 or 2.</div> </div>
Returned	<error code>[term]
Format	n
	<error code> Heater error code: 0 = no error, 1 = heater open load, 2 = heater short.
Remarks	Error condition is cleared upon querying the heater status, which will also clear the front panel error message.
<hr/>	
<b>IEEE</b>	<b>IEEE-488 Interface Parameter Command</b>
Input	<b>IEEE</b> <address>[term]
Format	nn
	<address> Specifies the IEEE address: 1–30. (Address 0 and 31 are reserved.)
Example	<b>IEEE 4[term]</b> after receipt of the current terminator, the instrument responds to address 4.

<b>IEEE?</b> Input Returned Format	<b>IEEE-488 Interface Parameter Query</b> <b>IEEE?</b> [term] <address>[term] nn (refer to command for description)
<b>INCRV</b> Input Format	<b>Input Curve Number Command</b> <b>INCRV</b> <input>,<curve number>[term] a,nn <input> Specifies which input to configure: A or B. <curve number> Specifies which curve the input uses. If specified curve type does not match the configured input type, the curve number defaults to 0. Valid entries: 0 = none, 1–20 = standard curves, 21–59 = user curves.
Remarks Example	Specifies the curve an input uses for temperature conversion. <b>INCRV A,23[term]</b> Input A uses User Curve 23 for temperature conversion.
<b>INCRV?</b> Input Format	<b>Input Curve Number Query</b> <b>INCRV?</b> <input>[term] a <input> Specifies which input to query: A or B.
Returned Format	<curve number>[term] nn (refer to command for description)
<b>INNAME</b> Input Format	<b>Sensor Input Name Command</b> <b>INNAME</b> <input>,<name>[term] a,s[15] <input> Specifies input to configure: A or B. <name> Specifies the name to associate with the sensor input.
Example	<b>INNAME A, "Sample Space"[term]</b> the string "Sample Space" will appear on the front panel display when possible to identify the sensor information being displayed.
Remarks	Be sure to use quotes when sending strings, otherwise characters such as spaces, and other non alpha-numeric characters, will be interpreted as a delimiter and the full string will not be accepted. It is not recommended to use commas or semi-colons in sensor input names as these characters are used as delimiters for query responses.
<b>INNAME?</b> Input Format	<b>Sensor Input Name Query</b> <b>INNAME?</b> <input>[term] a <input> Specifies input to query: A or B.
Returned Format	<name>[term] s[15] (refer to command for description)

<b>INTYPE</b>	<b>Input Type Parameter Command</b>	
Input	<b>INTYPE</b> <input>,<sensor type>,<autorange>,<range>,<compensation>,<units> [term]	
Format	a,n,n,n,n,n	
	<input>	Specifies input to configure: A or B
	<sensor type>	Specifies input sensor type: 0 = Disabled 1 = Diode 2 = Platinum RTD 3 = NTC RTD 4 = Thermocouple
	<autorange>	Specifies autoranging: 0 = off and 1 = on.
	<range>	Specifies input range when autorange is off:

Diode	0 = 2.5 V
	1 = 10 V
PTC RTD	0 = 10 $\Omega$
	1 = 30 $\Omega$
	2 = 100 $\Omega$
	3 = 300 $\Omega$
	4 = 1 k $\Omega$
	5 = 3 k $\Omega$
NTC RTD	6 = 10 k $\Omega$
	0 = 10 $\Omega$
	1 = 30 $\Omega$
	2 = 100 $\Omega$
	3 = 300 $\Omega$
	4 = 1 k $\Omega$
	5 = 3 k $\Omega$
	6 = 10 k $\Omega$
Thermocouple	7 = 30 k $\Omega$
	8 = 100 k $\Omega$
	0 = 50 mV

TABLE 6-8 Input range

	<compensation>	Specifies input compensation where 0 = off and 1 = on. Reversal for thermal EMF compensation if input is resistive, room compensation if input is thermocouple. Always 0 if input is a diode.
	<units>	Specifies the preferred units parameter for sensor readings and for the control setpoint: 1 = kelvin, 2 = Celsius, 3 = Sensor
Example	<b>INTYPE A,2,1,0,1,1[term]</b> sets Input A sensor type to Platinum RTD, autorange on, thermal compensation on, and preferred units to kelvin.	
Remarks	The <autorange> and <range> parameters do not apply to thermocouple sensor type, and the <autorange> and <compensation> parameters do not apply to diode sensor type. When configuring diode or thermocouple sensor types, these parameters must be included, but are ignored. A setting of 0 for each is recommended in this case.	

<b>INTYPE?</b>	<b>Input Type Parameter Query</b>	
Input	<b>INTYPE?</b> <input>[term]	
Format	a	
	<input>	Specifies input to query: A or B.
Returned	<sensor type>,<autorange>,<range>,<compensation>,<units> [term]	
Format	n,n,n,n,n (refer to command for description)	
Remarks	If autorange is on, the returned range parameter is the currently auto-selected range.	

<b>KRDG?</b>	<b>Kelvin Reading Query</b>
Input	<b>KRDG?</b> <input> [term]
Format	a <input> Specifies which input to query: A–B.
Returned	<kelvin value>[term]
Format	±nnnnnn
Remarks	Also see the RDGST? query.
<b>LEDS</b>	<b>Front Panel LEDS Command</b>
Input	<b>LEDS</b> <off/on> [term]
Format	n <off/on> 0 = LEDs Off, 1 = LEDs On
Remarks	If set to 0, front panel LEDs will not be functional. Function can be used when display brightness is a problem.
Example	<b>LED 0</b> [term] turns all front panel LED functionality off.
<b>LEDS?</b>	<b>Front Panel LEDS Query</b>
Input	<b>LEDS?</b> [term]
Returned	<off/on> [term]
Format	n (refer to command for description)
<b>LOCK</b>	<b>Front Panel Keyboard Lock Command</b>
Input	<b>LOCK</b> <state>, <code> [term]
Format	n,nnn <state> 0 = Unlocked, 1 = Locked <code> Specifies lock-out code. Valid entries are 000–999.
Remarks	Locks out all front panel entries except pressing <b>ALL OFF</b> to immediately turn off all heater outputs. Refer to section 4.7.
Example	<b>LOCK 1,123</b> [term] enables keypad lock and sets the code to 123.
<b>LOCK?</b>	<b>Front Panel Keyboard Lock Query</b>
Input	<b>LOCK?</b> [term]
Returned	<state>,<code>[term]
Format	n,nnn (refer to command for description)
<b>MDAT?</b>	<b>Minimum/Maximum Data Query</b>
Input	<b>MDAT?</b> <input> [term]
Format	a <input> Specifies which input to query: A or B.
Returned	<min value>,<max value>[term]
Format	±nnnnnn,±nnnnnn
Remarks	Returns the minimum and maximum input data. Also see the RDGST? query.
<b>MNMXRST</b>	<b>Minimum and Maximum Function Reset Command</b>
Input	<b>MNMXRST</b> [term]
Remarks	Resets the minimum and maximum data for all inputs.





<b>OPSTR?</b> Input Returned Format Remarks	<b>Operational Status Register Query</b> <b>OPSTR?</b> [term] <bit weighting> [term] nnn The integers returned represent the sum of the bit weighting of the operational status bits. These status bits are latched when the condition is detected. This register is cleared when it is read. Refer to section 6.2.5.2 for a list of operational status bits.
<b>OUTMODE</b> Input Format	<b>Output Mode Command</b> <b>OUTMODE</b> <output>, <mode>, <input>, <powerup enable> [term] n,n,n,n <output> Specifies which output to configure: 1 or 2. <mode> Specifies the control mode. Valid entries: 0 = Off, 1 = Closed Loop PID, 2 = Zone, 3 = Open Loop, 4 = Monitor out, 5 = Warmup Supply <input> Specifies which input to use for control: 0 = None, 1 = A, 2 = B <powerup enable> Specifies whether the output remains on or shuts off after power cycle. Valid entries: 0 = powerup enable off, 1 = powerup enable on.  <b>Example</b> <b>OUTMODE 1,2,1,0</b> [term] Output 1 configured for Zone control mode, using Input A for the control input sensor, and will turn the output off when power is cycled. <b>Remarks</b> Modes 4 and 5 are only valid for Output 2 in Voltage mode.
<b>OUTMODE?</b> Input Format Returned Format	<b>Output Mode Query</b> <b>OUTMODE?</b> <output> [term] n <output> Specifies which output to query: 1 or 2. <mode>, <input>, <powerup enable> [term] n,n,n (refer to command for description)
<b>PID</b> Input Format	<b>Control Loop PID Values Command</b> <b>PID</b> <output>, <P value>, <I value>, <D value> [term] n,+nnnnn,+nnnnn,+nnnn <output> Specifies which output's control loop to configure: 1 or 2. <P value> The value for output Proportional (gain): 0.1 to 1000. <I value> The value for output Integral (reset): 0.1 to 1000. <D value> The value for output Derivative (rate): 0 to 200.  <b>Remarks</b> Control settings, (P, I, D, and Setpoint) are assigned to outputs, which results in the settings being applied to any loop formed by the output and its control input. <b>Example</b> <b>PID 1,10,50,0</b> [term] Output 1 P is 10, I is 50, and D is 0%.
<b>PID?</b> Input Format Returned Format	<b>Control Loop PID Values Query</b> <b>PID?</b> <output> [term] n <output> Specifies which output's control loop to query: 1 or 2. <P value>, <I value>, <D value> [term] +nnnnn,+nnnnn,+nnnn (refer to command for description)

<hr/>	
<b>POLARITY</b>	<b>Output Voltage Polarity Command</b>
Input	<b>POLARITY</b> <output>,<polarity>[term]
Format	n,n <output>                      Output for which to configure the polarity setting: 2 <polarity>                  Specifies output voltage is 0=unipolar (positive output only) or 1=bipolar (positive or negative output)
Remarks	The polarity command only applies to Output 2, and only when output type is voltage.
<hr/>	
<b>POLARITY?</b>	<b>Output Voltage Polarity Query</b>
Input	<b>POLARITY?</b> [term]
Returned	<polarity>[term]
Format	n (refer to command for description)
<hr/>	
<b>RAMP</b>	<b>Control Setpoint Ramp Parameter Command</b>
Input	<b>RAMP</b> <output>,<off/on>,<rate value>[term]
Format	n,n,nnnn <output>                      Specifies which output's control loop to configure: 1 or 2. <off/on>                      Specifies whether ramping is 0 = Off or 1 = On. <rate value>                Specifies setpoint ramp rate in kelvin per minute from 0.1 to 100. The rate is always positive, but will respond to ramps up or down. A rate of 0 is interpreted as infinite, and will therefore respond as if setpoint ramping were off.
Example	<b>RAMP 1,1,10.5</b> [term]    when Output 1 setpoint is changed, ramp the current setpoint to the target setpoint at 10.5 K/minute.
Remarks	Control loop settings are assigned to outputs, which results in the settings being applied to the control loop formed by the output and its control input.
<hr/>	
<b>RAMP?</b>	<b>Control Setpoint Ramp Parameter Query</b>
Input	<b>RAMP?</b> <output>[term]
Format	n <output>                      Specifies which output's control loop to query: 1 or 2.
Returned	<off/on>,<rate value>[term]
Format	n,nnnn (refer to command for description)
<hr/>	
<b>RAMPST?</b>	<b>Control Setpoint Ramp Status Query</b>
Input	<b>RAMPST?</b> <output>[term]
Format	n <output>                      Specifies which output's control loop to query: 1 or 2.
Returned	<ramp status>[term]
Format	n <ramp status>                0 = Not ramping, 1 = Setpoint is ramping.

---

**RANGE Heater Range Command**

Input	<b>RANGE</b> <output>,<range>[term]
Format	n,n <output>                    Specifies which output to configure: 1 or 2. <range>                    For Outputs 1 and 2 in Current mode: 0 = Off, 1 = Low, 2 = Medium, 3 = High For Output 2 in Voltage mode: 0 = Off, 1 = On
Remarks	The range setting has no effect if an output is in the Off mode, and does not apply to an output in Monitor Out mode. An output in Monitor Out mode is always on.

---

**RANGE? Heater Range Query**

Input	<b>RANGE?</b> <output>[term]
Format	n <output>                    Specifies which output to query: 1 or 2.
Returned	<range>[term]
Format	n (refer to command for description)

---

**RDGST? Input Reading Status Query**

Input	<b>RDGST?</b> <input>[term]
Format	a <input>                    Specifies which input to query: A or B.
Returned	<status bit weighting>[term]
Format	nnn
Remarks	The integer returned represents the sum of the bit weighting of the input status flag bits. A "000" response indicates a valid reading is present.

Bit	Bit Weighting	Status Indicator
0	1	invalid reading
4	16	temp underrange
5	32	temp overrange
6	64	sensor units zero
7	128	sensor units overrange

---

**RELAY Relay Control Parameter Command**

Input	<b>RELAY</b> <relay number>,<mode>,<input alarm>,<alarm type>[term]
Format	n,n,a,n <relay number>            Specifies which relay to configure: 1 or 2. <mode>                    Specifies relay mode. 0 = Off, 1 = On, 2 = Alarms. <input alarm>            Specifies which input alarm activates the relay when the relay is in alarm mode: A or B. <alarm type>            Specifies the input alarm type that activates the relay when the relay is in alarm mode. 0 = Low alarm, 1 = High Alarm, 2 = Both Alarms.
Example	<b>RELAY 1,2,B,0</b> [term] relay 1 activates when Input B low alarm activates.

---

**RELAY? Relay Control Parameter Query**

Input	<b>RELAY?</b> <relay number>[term]
Format	n <relay number>            Specifies which relay to query: 1 or 2.
Returned	<mode>,<input alarm>,<alarm type>[term]
Format	n,a,n (refer to command for description)

<hr/>	
<b>RELAYST?</b>	<b>Relay Status Query</b>
Input	<b>RELAYST? &lt;relay number&gt;[term]</b>
Format	n <relay number>                      Specifies which relay to query: 1 or 2.
Returned	<status>[term]
Format	n                                      0 = Off, 1 = On.
<hr/>	
<b>SCAL</b>	<b>Generate SoftCal Curve Command</b>
Input	<b>SCAL &lt;std&gt;,&lt;dest&gt;,&lt;SN&gt;,&lt;T1 value&gt;,&lt;U1 value&gt;,&lt;T2 value&gt;,&lt;U2 value&gt;,&lt;T3 value&gt;,&lt;U3 value&gt;[term]</b>
Format	n,nn,S[10],+nnnnnn,±nnnnnn,±nnnnnn,±nnnnnn,±nnnnnn,±nnnnnn <std>                              Specifies the standard curve from which to generate a SoftCal™ curve. Valid entries: 1, 6, 7. <dest>                            Specifies the user curve to store the SoftCal™ curve. Valid entries: 21–59. <SN>                              Specifies the curve serial number. Limited to ten characters. <T1 value>                       Specifies first temperature point in kelvin. <U1 value>                       Specifies first sensor units point. <T2 value>                       Specifies second temperature point in kelvin. <U2 value>                       Specifies second sensor units point. <T3 value>                       Specifies third temperature point in kelvin. <U3 value>                       Specifies third sensor units point.
Remarks	Generates a SoftCal™ curve. Refer to Paragraph 5.3.
Example	<b>SCAL 1,21,1234567890,4.2,1.6260,77.32,1.0205,300.0,0.5189[term]</b> generates a three-point SoftCal™ curve from standard curve 1 and saves it in user curve 21.
<hr/>	
<b>SETP</b>	<b>Control Setpoint Command</b>
Input	<b>SETP &lt;output&gt;,&lt;value&gt;[term]</b>
Format	n,±nnnnnn <output>                           Specifies which output's control loop to configure: 1 or 2. <value>                           The value for the setpoint (in the preferred units of the control loop sensor).
Example	<b>SETP 1,122.5[term]</b> Output 1 setpoint is now 122.5 (based on its units).
Remarks	Control settings, that is, P, I, D, and Setpoint, are assigned to outputs, which results in the settings being applied to the control loop formed by the output and its control input.
<hr/>	
<b>SETP?</b>	<b>Control Setpoint Query</b>
Input	<b>SETP? &lt;output&gt;[term]</b>
Format	n <output>                           Specifies which output to query: 1 or 2.
Returned	<value>[term]
Format	±nnnnnn (refer to command for description)

<b>SRDG?</b>	<b>Sensor Units Input Reading Query</b>
Input	<b>SRDG?</b> <input> [term]
Format	a <input> Specifies which input to query: A or B.
Returned	<sensor units value>[term]
Format	±nnnnnn
Remarks	Also see the RDGST? command.
<b>TEMP?</b>	<b>Thermocouple Junction Temperature Query</b>
Input	<b>TEMP?</b> [term]
Returned	<junction temperature>[term]
Format	+nnnnn
Remarks	Temperature is in kelvin. This query returns the temperature of the ceramic thermocouple block used in the room temperature compensation calculation.
<b>TLIMIT</b>	<b>Temperature Limit Command</b>
Input	<b>TLIMIT</b> <input>, <limit> [term]
Format	n+nnnn <input> Specifies which input to configure: A or B. <limit> The temperature limit in kelvin for which to shut down all control outputs when exceeded. A temperature limit of zero turns the temperature limit feature off for the given sensor input.
Example	<b>TLIMIT B,450</b> [term] if the temperature of the sensor on Input B exceeds 450 K, all control outputs will be turned off.
Remarks	A temperature limit setting of 0 K turns the temperature limit feature off.
<b>TLIMIT?</b>	<b>Temperature Limit Query</b>
Input	<b>TLIMIT?</b> <input> [term]
Format	a <input> Specifies which input to query: A or B.
Returned	< limit>[term]
Format	+nnnn (refer to command for description)
<b>TUNEST?</b>	<b>Control Tuning Status Query</b>
Input	<b>TUNEST?</b> [term]
Returned	<tuning status>,<output>,<error status>,<stage status>[term]
Format	n,n,n,nn <tuning status> 0 = no active tuning, 1 = active tuning. <output> Heater output of the control loop being tuned (if tuning): 1 = Output 1, 2 = Output 2 <error status> 0 = no tuning error, 1 = tuning error <stage status> Specifies the current stage in the Autotune process. If tuning error occurred, stage status represents stage that failed.
Remarks	If initial conditions are not met when starting the autotune procedure, causing the autotuning process to never actually begin, then the error status will be set to 1 and the stage status will be stage 00.

**WARMUP Warmup Supply Parameter Command**

Input	<b>WARMUP</b> <output>,<control>,<percentage> [term]
Format	n,n,+nnn.nn <output>                      Output 2 is the only valid entry and must be included. <control>                    Specifies the type of control used: 0 = Auto Off, 1 = Continuous <percentage>                Specifies the percentage of full scale (10 V) Monitor Out voltage to apply to turn on the external power supply.
Example	<b>WARMUP 1,50[term]</b> Output 2 in Voltage mode will use continuous control, with a 5 V (50.50%) output voltage for activating the external power supply.
Remarks	Warmup mode applies only to Output 2 in Voltage mode. The Output Type parameter must be configured using the HTRSET command, and the Output mode and Control Input parameters must be configured using the OUTMODE command.

**WARMUP? Warmup Supply Parameter Query**

Input	<b>WARMUP?</b> <output> [term]
Returned	<control>,<percentage>[term]
Format	n,+nnn.nn (refer to command for description)

**ZONE Control Loop Zone Table Parameter Command**

Input	<b>ZONE</b> <output>,<zone>,<upper bound>,<P value>,<I value>,<D value>,<mout value>,<range>,<input>,<rate> [term]
Format	n,nn,+nnnnn,+nnnnn,+nnnnn,+nnnn,+nnnnn,n,n,+nnnn [term] <output>                      Specifies which heater output to configure: 1 or 2. <zone>                        Specifies which zone in the table to configure. Valid entries are: 1–10. <upper bound>                Specifies the upper Setpoint boundary of this zone in kelvin. <P value>                      Specifies the P for this zone: 0.1 to 1000. <I value>                      Specifies the I for this zone: 0.1 to 1000. <D value>                      Specifies the D for this zone: 0 to 200%. <mout value>                Specifies the manual output for this zone: 0 to 100%. <range>                        Specifies the heater range for this zone. Valid entries: 0 = Off, 1 = Low, 2 = Med, 3 = High. <input>                        Specifies the sensor input to use for this zone. 0 = Default (Use previously assigned sensor), 1 = Input A, 2 = Input B, <rate>                        Specifies the ramp rate for this zone: 0.1–100 K/min.
Remarks	Configures the output zone parameters. Refer to Paragraph 2.9.
Example	<b>ZONE 1,1,25.0,10,20,0,0,2,2,10[term]</b> Output 1 zone 1 is valid to 25.0 K with P = 10, I = 20, D = 0, a heater range of medium, sensor input B, and a ramp rate of 10 K/min.

**ZONE? Output Zone Table Parameter Query**

Input	<b>ZONE?</b> <output>,<zone> [term]
Format	n,nn <output>                      Specifies which heater output to query: 1 or 2. <zone>                        Specifies which zone in the table to query. Valid entries: 1–10.
Returned	< upper boundary>,<P value>,<I value>,<D value>,<mout value>,<range>,<input>,<rate>[term]
Format	+nnnnn,+nnnnn,+nnnnn,+nnnn,+nnnnn,n,n,+nnnn (refer to command for description)